

Digital Image Capture

The MediCap line of image capture devices allows you to easily save still images from any medical device that has standard video output. The USB200 model also allows you to save full motion video. Images and video are saved to a convenient USB flash drive.



MediCap USB200 - Video and Image Capture

The MediCap USB200 is the most powerful versatile solution for all modalities. This device easily captures full motion video or still images from your existing medical imaging equipments. It works with ultrasounds, endoscopes, intraoral cameras, surgical microscopes and most video



MediCap USB150 - Still Image Capture

The MediCap USB150 is an excellent solution for all modalities where only still images are needed. It easily captures digital images from medical imaging equipment such as ultrasounds, endoscopes, intraoral cameras, surgical microscopes and most video.



MediCap USB100 - Still Image Capture

The MediCap USB100 is an affordable solution for most modalities where only still images are needed. It easily captures digital images from medical imaging equipment such as ultrasounds, endoscopes, intraoral cameras, surgical microscopes and most video imaging devices.

Features and Benefits

- **Captures and stores digital images on USB flash drive.** Works with your existing medical imaging equipment.
- **Easy installation.** Just plug it in and capture.
- **Little or no training needed.** Simply push a button to capture the current image on your monitor - it is automatically saved to a USB flash drive.
- **Affordable.** Revolutionary technology allows the MediCap models to be priced at a level any practice, clinic or hospital can afford.
- **Small, portable and rugged.** The compact design easily fits on medical equipment carts and is designed for demanding hospital environments.
- **Powerful, user-friendly features.** Capture via button, foot switch or the trigger on your imaging device. Review your captured images on your monitor.

Comprehensive Solutions for Ultrasound Equipment

www.echoserve.com

Toll Free 877-386-3246

e
c
h
o
s
e
r
v
e

Digital Image Capture

MediCap Model Comparison

Features	USB100	USB150	USB200
Still Image Capture	✓	✓	✓
Video Capture			✓
USB Flash Drive Storage	✓	✓	✓
Included Flash Drive	256 MB	256 MB	512 MB
Remote Control / Foot Switch	✓	✓	✓
Built-In LCD Monitor		✓	✓
Built-In Speaker			✓
BNC Video Connection	✓	✓	✓
S-Video Connection	✓	✓	✓
800x600 Still Images	✓	✓	✓
1024x768 Still Images			✓
Conversion Required	✓	✓	
DICOM Still Images	✓	✓	✓
JPEG Still Images	✓	✓	✓
TIFF & PNG Still Images			✓
BMP & PCX Still Images	✓	✓	
MPEG-2 Video Files			✓
MPEG-1 Audio			✓
List Price	\$1,500	\$2,000	\$3,500

echoserve

Comprehensive Solutions for Ultrasound Equipment

www.echoserve.com

Toll Free 877-386-3246

Echoserve is an authorized distributor for MediCapture. Product names and images are used with permission of MediCapture.