

Electrical Leakage Testers



The ULT-2000 Series is specifically designed to test the electrical safety of all types of diagnostic ultrasound transducers independent of the ultrasound machine.

This specialty leakage tester has exact values for each manufacturer and probe type programmed into the unit, providing both pass/fail and quantitative values for each transducer.

Features Include:

- **Specific values for each probe model**
- Simple operation for non-technical users
- PC interface for data storage
- Software upgradeable for new probe types



The SA-2000 series is a microprocessor-based electrical safety analyzer. The primary function of these devices is to provide full electrical leakage testing for all types of biomedical equipment, including chassis leakage, transducer leakage, and ECG leakage tests.

Features Include:

- Touch control keys
- External resistance
- Externally replaceable ground fuse
- Built-in ECG simulator

All leakage testers feature compatibility with Dale Technology® DALE 800, Fluke Biomedical® ULT 800 adapters and electrodes allowing you to upgrade to the very latest technology platform, while safeguarding your prior investment in adapters and probes.

Comprehensive Solutions for Ultrasound Equipment

www.echoserve.com

Toll Free 1-877-386-3246

e
c
h
o
s
e
r
v
e

Leakage Tester Pricing

PN	Description	Price
ULT-2010	Ultrasound Probe Leakage Tester	\$950
ULT-2020	Ultrasound Probe Leakage Tester with Real-Time Monitoring and Logging	\$1,095
SA-2005	Electrical Safety Analyzer with 5-Lead ECG	\$1,080
SA-2010	Electrical Safety Analyzer with 10-Lead ECG	\$1,285
SA-2010S	Electrical Safety Analyzer with 10-Lead ECG and ECG Simulator	\$1,495
ULT-PC-10	Dual Conductivity Probe (Short Electrodes)	\$310
ULT-PC-20	Dual Conductivity Probe (Long Electrodes)	\$450
BC20-42200	CIDEX-Compatible Soak Tray	\$295
BC20-42300	External Printer	\$595
ULT-PA-11	Sonosite Probe Adapter	\$453
ULT-PA-10	Siemens Micropinless Probe Adapter	\$453
ULT-PA-12	Siemens 260-Pin Probe Adapter	\$515
ULT-PA-13	Siemens 156-Pin Probe Adapter	\$453
ULT-PA-14	Philips HDI 3000/3500/5000 Probe Adapter	\$448
ULT-PA-15	GE Logiq 7, Vivid 7, and Similar Probe Adapter	\$513
ULT-PA-16	GE Logiq 200 and Similar Probe Adapter	\$579
ULT-PA-17	Philips iE33/iU22 Probe Adapter	\$656
ULT-PA-18	Philips SONOS 4500/5500/7500 Probe Adapter	\$446
ULT-PA-20	Toshiba Probe Adapter	\$638
ULT-PA-21	Hitachi Probe Adapter	\$950
ULT-PA-22	Aloka Probe Adapter	\$525

ABOUT ECHOSERVE

Echoserve, Inc. is a private company headquartered in Golden, Colorado. We provide new and refurbished equipment, field service, technical support, training, and part depot repair/exchange services to hospitals, clinics, independent service organizations, and equipment dealers in the diagnostic ultrasound, general x-ray, and mammography markets. Our asset management solutions help customers minimize their total cost of ownership, maximize system up-time and achieve high levels of patient satisfaction.



1150 Catamount Drive, Golden, CO 80403 • 877-386-3246 • www.echoserve.com

Echoserve is an authorized distributor for BC Medical Group, International. Product names and images are used with permission of BC Medical Group.